


# Chamlide CMS

35mm dish type magnetic chambers for square coverslip


## 1-well Chamlide CMS




### Coverslip available

18 mm x 18 mm  
20 mm x 20 mm  
22 mm x 22 mm  
24 mm x 24 mm

- 1-well Chamlide CMS is for the square coverslip, and available on 18 x 18 / 20 x 20 / 22 x 22 / 24x24mm coverslip
- Other functions and features are the same as the standard type Chamlide CMB.




## 2-well Chamlide CMS




### Coverslip available

18 mm x 18 mm  
20 mm x 20 mm  
22 mm x 22 mm  
24 mm x 24 mm

- 2-well Chamlide CMS is for the square coverslip, and available on 18 x 18 / 20 x 20 / 22 x 22 / 24x24mm coverslip
- Other functions and features are the same as the standard type Chamlide CMM.




## 4-well Chamlide CMS



### Coverslip available

18 mm x 18 mm  
20 mm x 20 mm  
22 mm x 22 mm  
24 mm x 24 mm

- 4-well Chamlide CMS is for the square coverslip, and available on 18 x 18 / 20 x 20 / 22 x 22 / 24x24mm coverslip
- Other functions and features are the same as the standard type Chamlide CMM.



## Specifications

24x24mm coverslip type 1-well magnetic chamber internal volume

24x24mm coverslip type 2-well magnetic chamber internal volume

24x24mm coverslip type 4-well magnetic chamber internal volume

Chamber material

Maximum 2.0 ml / well

Maximum 800  $\mu\text{l}$  / well

Maximum 300  $\mu\text{l}$  / well

Body	Black PEEK
Bottom	Aluminum alloy
Cover	Glass
Gasket	Silicone

Model No.	Product
CM-S18-1	1 well Chamlide CMS for 18x18mm coverslip
CM-S20-1	1 well Chamlide CMS for 20x20mm coverslip
CM-S22-1	1 well Chamlide CMS for 22x22mm coverslip
CM-S24-1	1well Chamlide CMS for 24x24mm coverslip
CM-S18-2	2 well Chamlide CMS for 18x18mm coverslip
CM-S20-2	2 well Chamlide CMS for 20x20mm coverslip
CM-S22-2	2 well Chamlide CMS for 22x22mm coverslip
CM-S24-2	2well Chamlide CMS for 24x24mm coverslip
CM-S18-4	4 well Chamlide CMS for 18x18mm coverslip
CM-S20-4	4 well Chamlide CMS for 20x20mm coverslip
CM-S22-4	4 well Chamlide CMS for 22x22mm coverslip
CM-S24-4	4well Chamlide CMS for 24x24mm coverslip