


Well-Slip

Coverslip bottom magnetic chamber for rectangular coverslip




- Coverslip bottom multi-well chamber
- Replaces an expensive consumable chambered coverglass.
- Semi-permanent (re-useable).
- Simply remove and replace the bottom coverslip for next use.
- Autoclavable.
- Uses silicone gaskets to prevent leaking.
- Metal bottom for stable temperature control during live cell imaging.
- Non-toxic, non-reactive PEEK material.
- Patented magnetic system for easy assembly.
- Control the temperature and the CO₂ level using Chamlide TC for live cell imaging.



Specifications

Material	Main body	PEEK
	Bottom	Black anodized aluminum alloy
	Gasket	Silicone

Coverslip available	Quantities of wells
22 mm x 32 mm	1 well / 2 well
22 mm x 40 mm	1 well / 2 well / 4 well
22 mm x 50 mm	1 well / 2 well / 4 well / 8 well
22 mm x 60 mm	1 well / 2 well / 5 well / 10 well
24 mm x 32 mm	1 well / 2 well
24 mm x 40 mm	1 well / 2 well / 4 well
24 mm x 50 mm	1 well / 2 well / 4 well / 8 well
24 mm x 60 mm	1 well / 2 well / 5 well / 10 well



Well-slip for 22(24) x 60mm

Model No.	Product
WL-2232-1	1 Well-slip for 22x32mm coverslip
WL-2240-1	1 Well-slip for 22x40mm coverslip
WL-2250-1	1 Well-slip for 22x50mm coverslip
WL-2260-1	1 Well-slip for 22x60mm coverslip
WL-2432-1	1 Well-slip for 24x32mm coverslip
WL-2440-1	1 Well-slip for 24x40mm coverslip
WL-2450-1	1 Well-slip for 24x50mm coverslip
WL-2460-1	1 Well-slip for 24x60mm coverslip
WL-2232-2	2 Well-slip for 22x32mm coverslip
WL-2240-2	2 Well-slip for 22x40mm coverslip
WL-2250-2	2 Well-slip for 22x50mm coverslip
WL-2260-2	2 Well-slip for 22x60mm coverslip
WL-2432-2	2 Well-slip for 24x32mm coverslip

Model No.	Product
WL-2440-2	2 Well-slip for 24x40mm coverslip
WL-2450-2	2 Well-slip for 24x50mm coverslip
WL-2460-2	2 Well-slip for 24x60mm coverslip
WL-2240-4	4 Well-slip for 22x40mm coverslip
WL-2250-4	4 Well-slip for 22x50mm coverslip
WL-2440-4	4 Well-slip for 24x40mm coverslip
WL-2450-4	4 Well-slip for 24x50mm coverslip
WL-2260-5	5 Well-slip for 22x60mm coverslip
WL-2460-5	5 Well-slip for 24x60mm coverslip
WL-2250-8	8 Well-slip for 22x50mm coverslip
WL-2450-8	8 Well-slip for 24x50mm coverslip
WL-2260-10	10 Well-slip for 22x60mm coverslip
WL-2460-10	10 Well-slip for 24x60mm coverslip